



Developing 2D and 3D Micropolar FEM Models for Porous GBR Meshes in Dentistry Applications 22-24 October 2024, Florence, Italy



Introduction Guided Bone Regeneration (GBR) Mesh

- Using a mechanical barrier, such as a membrane, to separate and protect the area of bone loss from the surrounding tissue.
- 1 Exposing the Bone Defect: A small incision to access the area of defect.
- 2 Bone Grafting: a cement that may contain bone, antimicrobial additives and stimulants that promotes new bone growth is placed under the membrane.
- 3 Preparing GBR Mesh: cut and formed by the surgeon
- 4 Placing and fixing the GBR Mesh

COMSO

- New bone can grow in the designated area.
- Enough stiffness to create and maintain a suitable space for the intended bone regeneration.
- Enough porosity to facilitate the diffusion of fluids, oxygen, nutrients, and bioactive substances for cell growth.






2/26

Introduction Equivalent Homogenised Models for Porous Structures

- When the geometry includes pores, there is an internal structure.
- By using equivalent model and defining effective properties instead of the real microstructure, we can decrease computational expense of modelling and discretisation.
- Non-classical continuum theories can consider the internal structure without direct modelling of details.
- In the present work, non-classical micropolar theory as the homogenised equivalent model is implemented for the porous structures in 2D and 3D.





Introduction Micropolar Theory

COMSOL

October 22-24


ERENCE

- Degrees of Freedom: U, Ø
- Kinematic Relations: $E_{ij} = U_{j,i} + e_{jik} \Phi_k$

 $K_{ij} = \Phi_{j,i}$

- Balance Equations: $\Sigma_{ij,i} + P_j = 0$ $M_{ij,i} + e_{ijk} \Sigma_{jk} + Q_k = 0$
- Constitutive Equations: $\Sigma_{ij} = A_{ijkl}E_{kl} + B_{ijkl}K_{kl}$ (Linear Elastic) $M_{ii} = B_{iikl}E_{kl} + D_{ijkl}K_{kl}$

- **Φ** Micro-rotation
- *E* Nonsymmetric Strain
- K Curvature
- ∑ Nonsymmetric Stress
- M Couple-stress
- P Body Force
- Q Body Couple



- New strain measure as Curvature.
- New stress measure as Couple-stress.


Methodology 2D Tetragonal Micropolar Model for Porous Plate

- General linear elastic constitutive equations in plane
- Ortho-tetragonal Material Symmetry

CONFERENCE 2024 FLORENCE

October 22-24

 Considering the linear elastic micropolar in 2D, for the geometries we study here, 5 independent material parameters are required to define the equivalent model.




Methodology Homogenization of 2D Model

COMSO

October 22-24


- The heterogeneous structure, is called Micro Model and is described by Cauchy continuum.
- Homogenised model is called Macro Model and is described by constitutive parameters of a micropolar continuum.
- Kinematic Map: Links the two levels of description
- Express the microscopic displacement field within the RVE as a function of the macroscopic strain measures at the material point on the macro-level.

$$u^{\text{hom}} = E_{11} x + E_{12}^{SYM} y - \frac{K_2}{2} y^2 - K_1 xy - \frac{10}{L^2} \Theta(y^3 - 3yx^2)$$
$$v^{\text{hom}} = E_{12}^{SYM} x + E_{22} y + \frac{K_1}{2} x^2 + K_2 xy + \frac{10}{L^2} \Theta(x^3 - 3xy^2)$$



Methodology Homogenization of 2D Model


- Once the kinematic map is determined, we proceed to discover the micropolar material parameters.
- We use the finite element method (FEM) to compute the response of the porous structure under different loadings.
 - Extract the total elastic strain energy W_{FEM} stored in the RVE from COMSOL.
 - Equating to the energy of an equivalent micropolar continuum $W_{Micropolar}$ from analytical expression.
- For each scenario, the micropolar material parameters are determined in such a way that the homogenized model retains the same amount of strain energy when exposed to the same loading conditions.



 $M_1K_1 + M_2K_2$



Methodology Finding Effective Micropolar Material Parameters for 2D Model






 Equivalent homogenized model for Functionally Graded (FG) porous structure is derived by considering the homogenization procedure developed for homogenous porosities.

A1111, A1122 A1111, A1122

D11

A1212,A1221





()

 \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc

 $\bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc$

A1212,A1221

....

- First, a parametric study is conducted to find the equivalent parameters of uniform porous plates with various pore sizes.
- In this parametric study, the pore density (the number of pores per unit length) is kept constant.


0 0 0 0 0

A1212, A1221

D11

A1111, A1122 A1111, A1122

D11



Pore Size (*l p*)



D11

 \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc

 \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc

0 0 0 0 0

 $\bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc$

A1212,A1221

 Pore sizes are changed to find the required equivalent parameters for each section of the FG porous structure.

0 0 0 0 0


0 0 0 0 0

A1212,A1221

D11

A1111, A1122 A1111, A1122

D11



Pore Size (*l p*)



D11

 \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc

 \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc \bigcirc

0 0 0 0 0

A1212, A1221

These data are then utilised to find curve-fitted functions for each material parameter with respect to the pore size.

0 0 0


A1111, A1122 A1111, A1122

D11

0 0 0 0 0

A1212, A1221

D11



Pore Size (*l p*)





D11

A1212, A1221

Methodology Homogenization of 3D Model

- A size-effect is predicted by micropolar theory for the torsion of circular cylinders.
- It is possible to determine the micropolar parameters by measuring the torsional rigidity vs size.
- Specimens at different length scales, while keeping unit cell size constant.








Methodology **Homogenization of 3D Model**

- Bio-inspired 3D microstructures such as gyroids and other Triply Periodic Minimal Surfaces (TPMS) can be employed for designing GBR bone scaffolds.
- Studying the impact of changes in microstructural architecture which reflects in the equivalent micropolar parameters.
- The unit cell design, porous structure and CAD model are developed using nTop software.
- Final 3D meshed part is imported in COMSOL for further FEM analyses.

COMSOL

October 22-24



CAD Model

FEM Tests



Methodology FEM Micropolar Model in COMSOL Using PDE Weak Form

Multiplying balance equations by the test functions, integrating over computational domain D
 Strong Form

Balance Equations

COMSO

October 22-24

$$\Sigma_{ij,i} + P_j = 0$$
$$M_{ij,i} - e_{ijk} \Sigma_{ik} + Q_j = 0$$

$$\Rightarrow \int_{D} \Sigma_{ij,i} v_{U_j} + \int_{D} P_j v_{U_j} = 0$$

$$\Rightarrow \int_{D} M_{ij,i} v_{\Phi_j} - \int_{D} e_{ijk} \Sigma_{ik} v_{\Phi_j} + \int_{D} Q_j v_{\Phi_j} = 0$$

Weak Form
Test Function for Φ

And by using the divergence theorem:

Based on the product rule of derivatives:

$$\int_{D} \Sigma_{ij,i} v_{U_j} = \int_{D} \left(\Sigma_{ij} v_{U_j} \right)_{,i} - \int_{D} \left(\Sigma_{ij} v_{U_{j,i}} \right)$$
$$\int_{D} M_{ij,i} v_{\Phi_j} = \int_{D} \left(M_{ij} v_{\Phi_j} \right)_{,i} - \int_{D} \left(M_{ij} v_{\Phi_{j,i}} \right)$$

$$\begin{aligned} & -\int_{D} \Sigma_{ij} v_{U_{j,i}} + \int_{B} \Sigma_{ij} v_{U_{j}} n_{i} + \int_{D} P_{j} v_{U_{j}} = 0 \\ & -\int_{D} M_{ij} v_{\Phi_{j,i}} + \int_{B} M_{ij} v_{\Phi_{j}} n_{i} - \int_{D} e_{ijk} \Sigma_{ik} v_{\Phi_{j}} + \int_{D} Q_{j} v_{\Phi_{j}} = 0 \end{aligned}$$

Methodology FEM Micropolar Model in COMSOL Using PDE Weak Form

- By using partial differential equation (PDE) modelling in COMSOL instead of traditional FE modelling, no user subroutine is required.
- Various complex geometries, B.C., and loadings can be applied in a user-friendly graphical interface.
- Visualisation of the results is convenient.

- (s11*test(ux)+s12*test(uy)+s21*test(vx)+s22*test(vy)
+s23*test(vz)+s31*test(wx)+s32*test(wy)+s33*test(wz))
+ p1*test(u)+p2*test(v)+p3*test(w)


- (mu11*test(phi1x)+mu12*test(phi1y)+mu13*test(phi1z)
+mu21*test(phi2x)+mu22*test(phi2y)+mu23*test(phi2z)
+mu31*test(phi3x)+mu32*test(phi3y) + mu33*test(phi3z))
-((s12-s21)*test(phi3)+(s31-s13)*test(phi2)
+(s23-s32)*test(phi1))+q1*test(phi1)+q2*test(phi2)

+q3*test(phi3)


COMSOL

October 22-24

2024 FLORENCE




Methodology FEM Micropolar Model in COMSOL / Benchmarks for 2D Model






Methodology **FEM Micropolar Model in COMSOL / Benchmarks for 3D Model**





Developing 2D and 3D Micropolar FEM Models for Porous GBR Meshes in Dentistry Applications


Parametric Study for Various Pore Patterns Using COMSOL






Developing 2D and 3D Micropolar FEM Models for Porous GBR Meshes in Dentistry Applications 19/26

Parametric Study for Various Pore Patterns Using COMSOL / Parametrisation



COMSOL CONFERENCE 2024 FLORENCE October 22–24 Developing 2D

Parametric Study for Various Pore Patterns/Materials Using COMSOL






Discussion

Using the Equivalent Mechanical Parameters for Designing GBR Mesh

- **Optimum Design for Porous GBR Meshes:** Nearest mechanical properties to the bone.
- GBR mesh is in contact with the cortical (compact) bone.
- Try to make the material parameters of the GBR mesh consistent with its adjacent bone.
- Experimental estimate of the micropolar parameters of compact bone.
- Finding a configuration with the material parameters close to those reported for compact bone in the literature.



| | Value | Unit |
|--------------------------|---------------|------|
| A ₁₁₁₁ | 12.00 ~ 43.43 | GPa |
| A ₁₁₂₂ | 4.00 | GPa |
| A ₁₂₁₂ | 21.10 ~ 36.77 | GPa |
| A ₁₂₂₁ | -13.05 ~ 2.67 | GPa |
| D ₁₁ | 3.24 | kN |



COMSO

October 22-24

Experimental study of micropolar and couple stress elasticity in compact bone in bending J. Yang & R. Lakes


Discussion

Using the Equivalent Mechanical Parameters for Designing GBR Mesh

- Circular pores, Titanium Alloy
 - → Parameters are beyond the required values.
- Rectangular pores, Titanium Alloy

 \rightarrow Pore sizes of 0.13 L - 0.15 L and porosity about 0.7, a good agreement for A1122, A1212, A1221 can be achieved.

| | Value | Unit |
|--------------------------|---------------|------|
| A ₁₁₁₁ | 12.00 ~ 43.43 | GPa |
| A ₁₁₂₂ | 4.00 | GPa |
| A ₁₂₁₂ | 21.10 ~ 36.77 | GPa |
| A ₁₂₂₁ | -13.05 ~ 2.67 | GPa |
| D ₁₁ | 3.24 | kN |






Discussion

COMSO

Functionally Graded (FG) Porous Design for GBR Mesh


- GBR meshes are fixed to the mandible bone using biocompatible screws
- At these fixing locations higher stiffness is required. Therefore, smaller pore sizes are more desirable there.
- Mimicking natural FG structure of the bone
- FG structure of type O is suggested
- Part near fixing areas → as near as possible to cortical (compact) bone to provide required loadbearing capacities.





Discussion 3D FG Design for Porous GBR Dental Scaffold


- Currently, metallic GBR sheets are used as a mechanical barrier and bone grafts are inserted into the space created by the barrier.
- After the healing period, it is necessary another surgery for removing the GBR mesh.
- To avoid a second surgery, a recent solution is to use biodegradable materials such as polylactic acid (PLA) instead of metallic meshes.
- However, suboptimal mechanical properties of biodegradable materials such as low stiffness and strength have limited their application.






Discussion 3D FG Design for Porous GBR Dental Scaffold

- With the developed framework, an innovative 3D design for GBR meshes is suggested.
- The microstructure evolves from the porous structure at the lower level to the compact structure at the top surface.
- The porous structure can host bio-active agents for stimulating bone regeneration
- While the upper compact surface should provide the required stiffness.
- By implementing our model, the distribution in both the internal structure and the material properties (for instance using bimodal material such as PLA+HA (Hydroxy appetite) or adding nano-reinforcements) can be customized.



3D Design for GBR Dental Scaffold







Developing 2D and 3D Micropolar FEM Models for Porous GBR Meshes in Dentistry Applications

Questions



Developing 2D and 3D Micropolar FEM Models for Porous GBR Meshes in Dentistry Applications

Effectiveness of Homogenized Micropolar Model

- Comparison of displacement magnitude contours for detailed and homogenized micropolar model.
- Indentation Test.
- Predicted values are close to those of the detailed model

